



# Antlia Pro

## User Manual



ICALTECH INNOVATIONS PVT LTD

---

*Bangalore*

# Table of Contents

|                                                 |    |
|-------------------------------------------------|----|
| Table of Contents.....                          | 1  |
| List of Figures.....                            | 3  |
| Abbreviations .....                             | 5  |
| 1. Introduction and Features.....               | 6  |
| 1.1. Overview and Purpose of Antlia Pro .....   | 6  |
| 1.1.1. Key Functional Features.....             | 6  |
| 1.1.2. Main Advantages .....                    | 7  |
| 1.2. Copyright and Legal information .....      | 7  |
| 2. Technical Specifications.....                | 10 |
| 2.1. Technical Specification – Antlia Pro ..... | 10 |
| 2.2. Device Label .....                         | 11 |
| 3. Hardware Platform.....                       | 12 |
| 3.1. Components Overview .....                  | 12 |
| 3.1.1. Device Side View .....                   | 12 |
| 3.1.2. Device Rear View .....                   | 15 |
| 3.1.3. List of Accessories .....                | 16 |
| 3.2. Measurement Principle .....                | 17 |
| 3.2.1. Resistance (R) .....                     | 18 |
| 3.2.2. Reactance (X) .....                      | 19 |
| 3.2.3. Outcome Parameters.....                  | 19 |
| 4. Environmental Conditions .....               | 22 |
| 4.1. Operating Conditions .....                 | 22 |
| 4.2. Storage Conditions .....                   | 22 |
| 4.3. Transport Conditions .....                 | 23 |
| 5. An Overview .....                            | 24 |
| 5.1. Device Overview .....                      | 24 |
| 5.2. Intended users.....                        | 24 |
| 5.3. Installation and Setup .....               | 25 |
| 5.3.1. Unpacking & Inspection.....              | 25 |
| 5.3.2. Parking Position .....                   | 25 |
| 5.3.3. Installation .....                       | 26 |
| 5.3.4. Positioning .....                        | 26 |

|           |                                         |           |
|-----------|-----------------------------------------|-----------|
| <b>6.</b> | <b>Operating Instructions .....</b>     | <b>28</b> |
| 6.1.      | Safety Information.....                 | 28        |
| 6.1.1.    | General Safety Instructions.....        | 28        |
| 6.1.2.    | Equipment Safety.....                   | 29        |
| 6.1.3.    | Warning and Caution Notes .....         | 29        |
| 6.1.4.    | Contraindications.....                  | 32        |
| 6.1.5.    | Symbols.....                            | 33        |
| 6.2.      | Operational Training.....               | 35        |
| 6.3.      | Operating Procedure.....                | 35        |
| 6.3.1.    | User Login and Device Connection .....  | 35        |
| 6.3.2.    | Device Verification Test .....          | 37        |
| 6.3.3.    | Performing Patient Tests.....           | 39        |
| 6.4.      | Care And Maintenance .....              | 55        |
| 6.5.      | Troubleshooting .....                   | 57        |
| <b>7.</b> | <b>Application Software Design.....</b> | <b>63</b> |
| 7.1.      | Device connection with APP .....        | 63        |
| 7.1.1.    | Software Installation .....             | 63        |
| 7.1.2.    | Software Upgradation.....               | 67        |
| 7.1.3.    | Firmware Upgradation.....               | 68        |
| <b>8.</b> | <b>Exit from Application .....</b>      | <b>71</b> |
| 8.1.      | Software Closure.....                   | 71        |
| 8.2.      | Antlia Pro Device Power Off.....        | 71        |
|           | Warranty Certificate .....              | 72        |

## List of Figures

|                                                                            |    |
|----------------------------------------------------------------------------|----|
| Figure 1: Antlia Pro Device .....                                          | 12 |
| Figure 2: Rear View of Antlia Pro Device .....                             | 13 |
| Figure 3: LED Indication Panel .....                                       | 14 |
| Figure 4: Adjustable arm of the device .....                               | 15 |
| Figure 5: Rear View of Antlia Pro .....                                    | 16 |
| Figure 6: Graphical display of respiratory volume and mouth pressure ..... | 18 |
| Figure 7: Resistance(R) - Graphical Representation.....                    | 19 |
| Figure 8: Reactance (X) - Graphical Representation .....                   | 19 |
| Figure 9: Transport Position of Antlia Pro .....                           | 25 |
| Figure 10: Positioning of the Device and Laptop .....                      | 26 |
| Figure 11: Connecting the power supply adaptor .....                       | 27 |
| Figure 12: Device in Ready to use mode .....                               | 36 |
| Figure 13: Antlia Pro Application Login Page.....                          | 36 |
| Figure 14: Antlia Pro Home Screen .....                                    | 37 |
| Figure 15: Device Connection Status on Home Screen .....                   | 37 |
| Figure 16 Reference Impedance connected to FOT tube .....                  | 38 |
| Figure 17: Device verification - Impedance Value Check .....               | 38 |
| Figure 18: Patient Information Page .....                                  | 39 |
| Figure 19: Patient Information Entry Page .....                            | 40 |
| Figure 20: Positioning of Antlia Pro for Patient Test .....                | 40 |
| Figure 21: iCaltech Filter Connection.....                                 | 41 |
| Figure 22: Patient Positioning.....                                        | 41 |
| Figure 23: Patient Database on Conduct Test Page.....                      | 43 |
| Figure 24: Selecting Patient for Testing.....                              | 43 |
| Figure 25: Conduct Test Page .....                                         | 44 |
| Figure 26: Frequency Mode Selection .....                                  | 44 |
| Figure 27: Start Breathing Test Instruction.....                           | 45 |
| Figure 28: Volume and Mouth Pressure Measurement .....                     | 45 |
| Figure 29: Test Report Page.....                                           | 46 |
| Figure 30: Real Time Graph of Volume and Mouth Pressure Measurement .....  | 46 |
| Figure 31: Report screen Page .....                                        | 47 |
| Figure 32: Saving Report Saving to the PC .....                            | 47 |
| Figure 33: Antlia Pro Installation File.....                               | 63 |
| Figure 34: PC's Properties.....                                            | 64 |
| Figure 35: System Type - Verification .....                                | 64 |
| Figure 36: Installation Folder Location in PC .....                        | 65 |
| Figure 37: Selecting for Desktop Shortcut Icon.....                        | 66 |
| Figure 38: Installation Process - Completed .....                          | 66 |
| Figure 39: Antlia Pro - Desktop Icon .....                                 | 67 |
| Figure 40: Antlia pro - Main application screen .....                      | 67 |
| Figure 41: Data base file location in program files.....                   | 68 |
| Figure 42: Firmware Upgradation Page.....                                  | 68 |

|                                                               |    |
|---------------------------------------------------------------|----|
| Figure 43: Firmware upgrade mode option.....                  | 69 |
| Figure 44: Choose Firmware HEX File .....                     | 70 |
| Figure 45: Start Upgrade option on firmware upgrade page..... | 70 |

## Abbreviations

|        |                               |
|--------|-------------------------------|
| Rrs    | Respiratory Resistance        |
| Xrs    | Respiratory Reactance         |
| R5     | Resistance at 5Hz             |
| R11    | Resistance at 11 Hz           |
| R19    | Resistance at 19Hz            |
| X5     | Reactance at 5Hz              |
| X11    | Reactance at 11 Hz            |
| X19    | Reactance at 19Hz             |
| Z5     | Respiratory impedance at 5Hz  |
| Z19    | Respiratory impedance at 19Hz |
| AX     | Area of Reactance             |
| Fres   | Resonant frequency (Hz)       |
| R5 COV | R5 Coefficient of Variation   |
| Z5 COV | Z5 Coefficient of Variation   |
| BVF    | Bacterial Viral Filter        |
| PFT    | Pulmonary Function Test       |
| NA     | Not Applicable                |

# 1. Introduction and Features

## 1.1. Overview and Purpose of Antlia Pro

the Antlia Pro – a cutting-edge medical device utilizing the Forced Oscillation Technique (FOT) to measure respiratory system impedance. Designed for pediatric and adult patients aged 4 years and above, the Antlia Pro is an invaluable tool in evaluating and managing respiratory conditions such as asthma and chronic obstructive pulmonary disease (COPD).

By precisely measuring airway resistance, reactance, and other lung function parameters, the device aids healthcare professionals in diagnosing and treating respiratory diseases effectively in both children and adults.

### 1.1.1. Key Functional Features

The key features of the Antlia Pro device are designed to ensure accurate measurement and ease of use for healthcare professionals:

- Captures resistance and reactance across a frequency range of 5–37 Hz.
- Operates in two frequency modes:
  - Mono-frequency mode (5 Hz and 19 Hz).
  - Multi-frequency mode (5–37 Hz).
- Performs pulmonary function tests (PFT):
  - 12 seconds in Mono-frequency mode.
  - 24 seconds in Multi-frequency mode.
- Supports both PRE and POST test types for bronchodilator or bronchoconstrictor reversibility testing.
- Conducts system self-checks to ensure proper functionality:
  - Pressure Sensor offset.
  - USB Communication port.
- Factory-calibrated for impedance measurements in Mono-frequency mode.

### 1.1.2. Main Advantages

- Measures lung function with high sensitivity, capturing small and peripheral airway diseases.
- Useful for evaluating bronchodilator reversibility and distinguishing normal from abnormal lung function.
- Suitable for patients unable to perform spirometry tests (e.g., pediatric, geriatric, or severely ill patients).
- Assists in diagnosing conditions such as COPD, asthma, interstitial lung disease (ILD), fibrosis, airway clearance issues, and silicosis.
- Efficiently identifies rapid changes in impedance during the respiratory cycle, particularly in Mono-frequency mode.

## 1.2. Copyright and Legal information

This document contains proprietary information protected by copyright. It is strictly prohibited to copy, reproduce, translate, or distribute this manual into another language without the express written permission of iCaltech Innovations Pvt Ltd.

All other names, company names, product names, and trademarks mentioned in this manual are the property of their respective owners.

iCaltech follows a policy of continuous development and innovation. iCaltech Innovations Pvt Ltd reserves the right to make changes and improvements to any product or component described in this document without prior notice.

### Model Information

- **Model Name:** Antlia Pro
- **Product Code:** ANTPRO
- **Part Number:** 205-0001-000



## Contact Information

For additional information or assistance, please contact the manufacturer.

**Manufacturer:** iCaltech Innovations Pvt Ltd.

**Address:** 2nd Floor, Sy no 80, Municipal No 4 & 4A (Katha no 351),

Muneshwara Layout, Service Road, Ring Road,

Laggere Bridge, Laggere,

Bengaluru, Karnataka, 560058, India.

**Email:** [operations@iCaltech.com](mailto:operations@iCaltech.com)

**Toll-Free Number:** +91 9480017633

## Scope of User Manual

The **Antlia Pro User Manual** provides essential guidance for the unpacking, installation, configuration, safe operation, maintenance, and troubleshooting of the device. It ensures safe and efficient usage to achieve optimal performance in clinical environments.

This manual is intended for healthcare professionals and trained personnel responsible for assessing and monitoring respiratory conditions in patients. By covering all necessary aspects of setup and usage, it enables users to maximize the device's capabilities.

## About this User Manual

- **Name of Document:** Antlia Pro User Manual
- **Revision:** 05
- **Part Number:** 203-0001-004

We strongly recommend that you read the entire user manual thoroughly before using the device.

On behalf of the iCaltech team, our partners, and representatives, we sincerely thank you for placing your trust in us. We are dedicated to supporting you in maximizing the potential of this state-of-the-art technology to deliver exceptional care to your patients.

**Important Notes**

- Any patient or equipment data used as examples in this User Manual are entirely fictional and intended for representation purposes only.
- The actual product may vary slightly in appearance or features, depending on the model, installed options, and specific configuration.
- Any resemblance to real persons, institutions, or names is purely coincidental.

All parameters, graphs, and images shown in this User Manual are for illustrative purposes only.

The parameters displayed by your device during operation are definitive.

## 2. Technical Specifications

### 2.1. Technical Specification – Antlia Pro



# Technical Specification – Antlia Pro

|                                                                  |                                                                                                        |                                         |                |
|------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-----------------------------------------|----------------|
| Flow Measurement                                                 | Calorimetric                                                                                           | Physical Dimensions                     | Height: 230 mm |
| Pressure range                                                   | ± 2.5 kPa                                                                                              |                                         | Width: 140 mm  |
| Flow Range                                                       | ± 4.16 L/s                                                                                             |                                         | Length: 240 mm |
| Pressure Transducer                                              | Piezo Resistive                                                                                        | Weight                                  | 2.90 Kgs       |
| Testing Signal Mode                                              | Mono Frequency                                                                                         | Mixed Frequency (5 Hz – 37 Hz)          |                |
| Frequency Range                                                  | 5 and 19 Hz                                                                                            | 5, 11, 13, 17, 19, 23, 29, 31 and 37 Hz |                |
| Measurement Accuracy                                             | For Impedance ±10%                                                                                     |                                         |                |
| Calibration                                                      | Factory Calibration and Auto-Zeroing of the sensors and A calibration check with a Reference Impedance |                                         |                |
| Connectivity                                                     | USB Port                                                                                               |                                         |                |
| Applied Port                                                     | Type B                                                                                                 |                                         |                |
| Portability                                                      | Table Top                                                                                              |                                         |                |
| Data Viewing and Storage                                         | Desktop/ Laptop                                                                                        |                                         |                |
| Support Arm                                                      | Metal                                                                                                  |                                         |                |
| Input Voltage                                                    | 100 – 240 V AC, 1.3 – 0.6 A, 50/60Hz                                                                   |                                         |                |
| Output Voltage                                                   | 12v DC, 800 milli Amps to 1 Amps                                                                       |                                         |                |
| Power Consumption                                                | 36 W                                                                                                   |                                         |                |
| Design Comply with IEC 60601-1, IEC 60601-1-2 and ERS Guidelines |                                                                                                        |                                         |                |

## 2.2. Device Label



### 3. Hardware Platform



Figure 1: Antlia Pro Device

#### 3.1. Components Overview

##### 3.1.1. Device Side View

Main components of the Antlia Pro Unit:

- A. FOT tube enclosure
- B. Source tube
- C. Source Enclosure
- D. LED indication Panel
- E. Adjustable Arm



Figure 2: Rear View of Antlia Pro Device

**A. FOT tube Enclosure:**

The FOT tube enclosure consists of three parts: Top cover, Bottom cover, and Front cover. Any component assembled in the FOT setup can be accessed by removing the top cover.

**B. Source Tube:**

Antlia Pro device has a 50 cm long source tube, which is used to carry the pressure wave from Source to FOT tube.

**C. Source Enclosure:**

Antlia Pro includes a Source that generates pressure signals ranging from 5 Hz to 37 Hz, and a control unit that measures volume and pressure.

**D. LED Indication Panel:**

The source enclosure also houses the LED indication slot, which contains three LEDs: Green, Blue and Red, to notify the user in different scenarios.





Figure 3: LED Indication Panel

### LED status


| Color | Status      | Notification                                                                                      |
|-------|-------------|---------------------------------------------------------------------------------------------------|
| UNLIT | OFF         | The device is powered down                                                                        |
| Green | Blinking    | The device is powered up                                                                          |
| Green | Solid state | The device is connected to App                                                                    |
| Blue  | Blinking    | The device is in start mode                                                                       |
| Blue  | Solid state | Conducting test                                                                                   |
| Red   | Error       | The device is powered up, but a critical error has occurred, and the system is not ready for use. |

**E. Adjustable Arm**

Using the adjustable arm, we can adjust the height and the inclination of the device to suit the patient's height, allowing to perform test easily. Move the arm in the upward and downward directions to adjust the height of the FOT section.



Min. height of the device - 20.3cm



Max. height of the device - 49.5 cm

Figure 4: Adjustable arm of the device

**3.1.2. Device Rear View**

The main components of the Antlia Pro Device located on the rear side of device are as follows:

**A. Power ON/OFF button**

Switch to turn the device power ON or OFF.

**B. Power supply connector**

Insert power adaptor in to this connector to supply power to the device.

**C. Type B USB port**

Facilitates connection between the Antlia Pro device and the computer system.






Figure 5: Rear View of Antlia Pro

### 3.1.3. List of Accessories

| Component Description                                                                                       | Image                                                                                 |
|-------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>iCaltech Filter:</b><br>The BVF filter recommended by iCaltech, is intended for single-patient use only. |   |
| <b>Nose clip:</b><br>To close the nose during the test                                                      |  |
| <b>Power Adapter:</b><br>Desktop AC Adapters 36 W, 12 V, 3A Medical.                                        |  |

|                                                                                                                                                        |                                                                                     |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <p><b>Reference impedance:</b></p> <p>It is used for device verification and provides a known impedance with a tolerance of <math>\pm 10\%</math>.</p> |  |
| <p><b>Device cap:</b></p> <p>It is used to cover the FOT inlet, preventing the ingress of dust and humidity into the device when not in use.</p>       |  |

**Note:** To interface the Antlia Pro device with a laptop or computer system, it is recommended to use a Type-B USB cable that is less than 1 meter in length.

## 3.2. Measurement Principle

The device generates pressure signals against which the patient performs tidal breathing. It provides respiratory parameters, such as resistance and reactance, at mixed and mono frequencies to assist in diagnosing respiratory disorders. The real time graphical display of respiratory volume and mouth pressure enables the user to observe and identify the subject's breathing maneuver during the test.

An electric signal is generated in the impedance analyzer and applied to a pump, which converts it into a corresponding pressure waveform. This pressure wave typically has an amplitude of approximately 0.2 kPa (peak-to-peak).




Figure 6: Graphical display of respiratory volume and mouth pressure

- This wave is directed to the patient's mouth, via iCaltech filter. The airflow ( $V'$ ) is measured using pneumotachograph, and the pressure ( $P$ ) is measured at the patient's airway input. The resulting signals are filtered and digitized in the impedance analyzer, which then calculates the respiratory impedance ( $Z_{rs}$ ).
- The mechanical behavior of the human respiratory system depends on the frequency of the applied oscillation. Low frequencies provide insights in to the mechanical properties of lung periphery, while high frequencies describe the properties of central and upper airways. The Antlia Pro utilizes oscillatory waveform ranging from 5 Hz to 37 Hz.

### 3.2.1. Resistance (R)

Resistance describes the dissipative mechanical properties of the airways, as well as the lung and chest wall tissue. Resistance mainly reflects the frictional pressure losses that occur as air moves in and out the airways, making it closely related to the overall airway's diameter.




Figure 7: Resistance(R) - Graphical Representation

### 3.2.2. Reactance (X)

Reactance is related to the energy storage capacity of the respiratory system. It is determined jointly by the elastic properties (the relationship between pressure and volume), which are dominant at low oscillation frequencies, and the inertive properties (the relationship between pressure and volume acceleration), which are dominant at high frequencies.




Figure 8: Reactance (X) - Graphical Representation

### 3.2.3. Outcome Parameters

#### 3.2.3.1. FOT Parameters

| Parameter       | Abbreviation |                                                                  |
|-----------------|--------------|------------------------------------------------------------------|
| Resistance 5Hz  | R5           | Low-frequency Resistance (5 Hz) reflects total airway resistance |
| Resistance 19Hz | R19          | Resistance (19 Hz) reflects mostly conducting airways            |

|                                         |           |                                                                                                  |
|-----------------------------------------|-----------|--------------------------------------------------------------------------------------------------|
| Frequency dep. of Resistance, 5 – 19 Hz | R5-R19    | Reflects Resistance in the small airways                                                         |
| Reactance 5Hz                           | X5        | Low-frequency Reactance (5Hz) reflects elastic compliance & small airways                        |
| Resonance frequency                     | Fres (Hz) | Resonant frequency reflects peripheral airways.                                                  |
| Delta X5Hz                              | Delta X5  | Low frequency Delta (5Hz) reflects change between inspiratory reactance and expiratory reactance |
| Area of Reactance                       | AX        | Area of reactance reflects elastic properties of lungs                                           |
| R5 COV                                  | R5 COV    | R5 Coefficient of variation                                                                      |
| Z5 COV                                  | Z5 COV    | Z5 Coefficient of variation                                                                      |

### 3.2.3.2. Breathing Parameters

| Parameter             | Abbreviation |                                                                             |
|-----------------------|--------------|-----------------------------------------------------------------------------|
| X5 Inspiration        | X5Insp       | Reactance value at inspiration cycle                                        |
| X5 Expiration         | X5Exp        | Reactance value at expiration cycle                                         |
| R5 Inspiration        | R5Insp       | Resistance value at inspiration cycle                                       |
| R5 Expiration         | R5Exp        | Resistance value at expiration cycle                                        |
| Tidal Volume          | Tidal Volume | Amount of air that moves in or out of the lungs with each respiratory cycle |
| Time Inspiration      | TimeInsp     | Average time duration for inspiration                                       |
| Time Expiration       | TimeExp      | Average time duration for expiration                                        |
| Respiratory Rate      | RR           | Number of breaths per minute                                                |
| R5 End Inspiration    | R5EndInsp    | Resistance at end of the inspiration cycle                                  |
| R5 End Expiration     | R5EndExp     | Resistance at end of the expiration cycle                                   |
| X5 End Inspiration    | X5EndInsp    | Reactance at end of inspiration cycle                                       |
| X5 End Expiration     | X5EndExp     | Reactance at end of expiration cycle                                        |
| Mean Inspiration Flow | Qi           | Average inspiratory flow rate                                               |

|                      |       |                                                                     |
|----------------------|-------|---------------------------------------------------------------------|
| Mean Expiration Flow | $Q_e$ | Average expiratory flow rate                                        |
| Minutes Ventilation  | $V_e$ | Volume of gas inhaled and exhaled from a person's lungs per minute. |

## 4. Environmental Conditions

The suitable environmental conditions for the safe operation, storage, and transport of the Antlia Pro are outlined below, ensuring accurate test results and promoting the long life of the product.

### 4.1. Operating Conditions

- **Air Quality:** The device should be used in a clean, dry environment free from dust and away from explosive and flammable substances.
- **Ventilation:** Ensure adequate space around the device to prevent overheating and for the patient accessibility. Keep the device away from the heat sources.
- **Electromagnetic Interference:** Do not operate the device near the strong electrical and magnetic devices such as X Ray Machines, Transformer, or Power lines to ensure accurate performance.
- **Vibration:** Do not use device near highly vibrating equipment like Generators, Large Motors, or in Rooms/Areas near to roads which has access to heavy vehicles.
- **Temperature:** +10°C to +40°C
- **Humidity:** 20% to 90% Relative Humidity (non-condensing)
- **Atmospheric Pressure:** 600 mbar (600 kPa) to 1060 mbar (106 kPa)
- **Altitude level:** Not to exceed 2,000 meters above sea level

### 4.2. Storage Conditions

- **Temperature:** 5°C to +40°C
- **Humidity:** 5% to 95% relative humidity (non-condensing)

### 4.3. Transport Conditions

- **Temperature:** 5°C to +40°C
- **Humidity:** 5% to 95% relative humidity (non-condensing)
- **Shock and Vibration:** Handle with care to avoid mechanical shock or excessive vibration during transit.
- **Protection from Exposure:** Avoid exposing the device to direct sunlight, rain, or extreme environmental conditions during transport.



## 5. An Overview

### 5.1. Device Overview

The Antlia Pro is a professional, PC-based device designed for the non-invasive measurement of the mechanical properties of the human respiratory system. It uses the Forced Oscillation Technique (FOT) to measure these properties.

The respiratory system with an oscillation pressure is used to evaluate the response in terms of mechanical impedance. Impedance represents the complex relationship between flow and mouth pressure.

Minimal patient cooperation is required during the test, as the device can capture the mechanical properties of the respiratory system while the patient breathes normally. It is suitable for both adult and pediatric patients (ages 4 years and above).

The test duration is brief: the mono-frequency test lasts 12 seconds, while the multi-frequency test lasts 24 seconds. During the test, the patient should breathe normally by holding their cheeks.

The Antlia Pro device should not operate continuously. Each "Act" (patient's breathing operation) lasts approximately 24 seconds, followed by a 2-minute rest period. This cycle ensures that patients do not experience discomfort due to the forced oscillations.

### 5.2. Intended users

The **Antlia Pro** is designed for use by healthcare professionals with expertise in respiratory care, including:

- **Pulmonologists**
- **Chest physicians**
- **General physicians**
- **Respiratory therapists**
- **Lab technicians**
- Other trained personnel in hospitals, clinics, and healthcare facilities.

The device is intended to assist professionals in accurately diagnosing, monitoring, and managing respiratory conditions, ensuring effective treatment and care for patients.

## 5.3. Installation and Setup

### 5.3.1. Unpacking & Inspection

- Before opening, inspect the Antlia Pro packaging for damage during shipping.
- If the box appears damaged, do not open the box and immediately contact iCaltech Technical Support or your local representative.
- Open the Antlia Pro device packaging carefully without allowing sharp objects to protrude into the box.
- Remove packing foam and carefully remove all items from the box.
- Visually inspect all items, and make sure that you obtained all system components as per Antlia Pro Package list.

### 5.3.2. Parking Position

When device is not in use or being transported, the device must be set to the transport position as shown in the below picture.



Figure 9: Transport Position of Antlia Pro

Note: During device transport, the Power Adapter cable, USB cable and Accessories shall be disconnected from the device.

### 5.3.3. Installation

Before starting the installation and operation of the device, make sure that the environmental condition reported in section 2.6 are met and that the laptop/ computer system configuration meets the following requirements:

- a) Processor: Intel i5 processor or above
- b) Operating System: Windows 10/11 - 64 bit or more
- c) RAM: 16 GB or more
- d) Free Space in Hard Drive: 2 GB or more
- e) Display: High Definition (HD) 1280 × 720

### 5.3.4. Positioning

- Place the Antlia Pro Unit on a desk or other flat surface within easy reach of the intended patient location during the measurement. Position the computer in such away, that the operator can easily view both the patient and the computer screen simultaneously.



Figure 10: Positioning of the Device and Laptop

- Connect the power supply adapter to its jack connector, which is located on the rear side of the device, and plug the adapter into the power socket.



Figure 11: Connecting the power supply adaptor

- Connect the PC to the unit using Type B - USB cable and Switch on the device by pressing power ON button, which is located on rear side of the device.
- In the power ON condition, the Green LED will blink on LED panel located on the Top panel of the unit, which is on the front side of the device.

## 6. Operating Instructions

Before operating the device, the operator should be aware of safety information about the device, including symbols marked on the product and package. They should also understand contraindications, warnings, and cautions that may arise during operation and their potential causes

### 6.1. Safety Information

#### 6.1.1. General Safety Instructions

To ensure the safety of patients, operators, and the device, follow these guidelines:

- Inspect power cords and plugs for damage before powering on the device. Replace damaged components through authorized personnel only.
- Keep this User Manual accessible for the users and adhere strictly to the instructions provided.
- Use only accessories and consumables recommended by iCaltech. Unapproved items may compromise safety and device performance.
- Protect all cables from heat, sharp objects, and rough surfaces. Regularly inspect cables for damage.
- Avoid exposing the device to direct sunlight or excessive light (over 1500 lux).
- Operate the device within the specified environmental conditions to ensure optimal performance.
- The device should not be used with other medical equipment unless explicitly approved by **iCaltech**.
- Ensure the breathing circuit is free from contaminants to prevent serious health risks, including aspiration and life-threatening conditions.
- Handle patient data as per the latest health data protection laws and regulations.

### 6.1.2. Equipment Safety

Improper handling of the Antlia Pro may lead to device damage, patient harm, or safety risks. **iCaltech** disclaims responsibility for:

- Operating the device in any manner other than as specified in this User Manual.
- Performing upgrades, adjustments, modifications, or repairs by personnel who are not contracted or authorized by iCaltech.
- Failing to replace components that affect the device's safe operation with original spare parts in the event of a malfunction.
- Using the device in an environment where the electrical wiring does not comply with local safety regulations.

### 6.1.3. Warning and Caution Notes

**Warning:** Indicates a potential hazard that may result in serious injury or death if ignored.

**Caution:** Indicates a hazard that may cause minor or moderate injury, or equipment damage if not observed.



#### **WARNING**

Do not open the device, as it contains no user serviceable parts inside.



#### **WARNING**

If any part of the device or its accessories such as Power adapter is changed, device performance may degrade performance and result in inaccurate results.



#### **WARNING**

##### **Improper use of the device:**

Do not expose the device to rain, moisture, or wet surfaces.



#### **WARNING**

If any part of the device or its accessories (e.g., power adapter, cables) is damaged or missing, discontinue use immediately. Using a damaged component could result in injury to the patient or operator.

**WARNING**

To mount the device, use a table of appropriate size and structure capable of supporting its weight.

**WARNING**

Unsuitable environmental conditions can affect the device's performance and lead to improper measurements. Use the device only under the specified environmental conditions.

**WARNING**

Do not modify the device in any way. Modifications can compromise safety and cause measurement errors.

**WARNING**

Do not position the device in a way that makes it difficult to disconnect the power cord in an emergency.

**WARNING**

The iCaltech filter and nose clip are intended for single use only.

**WARNING**

Do not perform service or maintenance activities while the device is in use with a patient.

**WARNING**

Do not use this device in the presence of flammable anesthetics.



### CAUTION

#### **Operation of the Device by Untrained or Unauthorized Users**

Risk of incorrect measurements due to the misinterpretation of system information.

The system should only be operated by individuals with proper training.



### CAUTION

#### **Use of Unapproved Software or Hardware Components**

Risk of malfunctions that could pose a hazard to patients and equipment.

Only use software or hardware components provided by iCaltech. The use of unauthorized accessories will void the warranty.



### CAUTION

#### **Modification to the System**

Risk of system malfunction.

Users are not authorized to change the calibration or modify/upgrade any software. Please contact iCaltech service team for assistance.



### CAUTION

#### **Pulling/Pinching the Source Tube**

Risk of incorrect measurements or results.

Users should not damage the source tube. The device should operate by trained persons, and care should be taken during transportation.

If the source tube is damaged, please contact iCaltech service team for assistance.



### CAUTION

#### **Blockage of Inertance Tube**

Risk of Malfunction.



Users must not block the inertance tube opening. Before conducting tests, ensure that there is no blockage near the opening of the inertance tube.

**CAUTION****Electromagnetic Disturbance**

If the unit hangs or resets due to electromagnetic disturbance, the operator conducting the test must reconnect the device to the app and perform a new test with the patient. Since multiple breathing tests are required from the patient, and tidal breathing (normal breathing) is needed, this degradation is acceptable.

**CAUTION****System Date and Time**

The system date and time on the computer must be accurate. If not, please update it to ensure that the date and time on the report are printed correctly, as they are derived from the computer's settings.

**CAUTION**

Application stops the test if the device runs continuously in start mode for 60 seconds.

The Test stops, if no patient breathing is detected.

**CAUTION**

Use only iCaltech recommended filters. Using any other filters may lead to incorrect results and impact the safe operation of the product.

### 6.1.4. Contraindications

Device is not supposed to be used for patients with any of the following condition:


- Acute disorders affecting test performance (e.g., Vomiting, nausea, vertigo)
- Recent eye (for example, cataract), thoracic and abdominal surgery
- Oral or facial pain exacerbated by a BV Filter.


- Recent myocardial infarction

Additionally, please do not operate the device.

- In CT and MRI environment.
- In OT environment.
- Along with high surgical equipment.

### 6.1.5. Symbols

| Symbol                                                                              | Description/ Meaning of Symbol                                                                    |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|
|    | This symbol indicates that the device operates on DC current.                                     |
|    | General warning sign                                                                              |
|    | This symbol instructs the user to consult the device's instructions for proper usage information. |
|   | Accompanied by a date, this symbol indicates the device's manufacturing date.                     |
|  | Serial number. This symbol should be accompanied by the device's serial number.                   |
|  | Refers the legal Manufacturer (iCaltech Innovations Pvt Ltd)                                      |
|  | PE B applied part: iCaltech Filter and Nose clip are applied parts                                |
|  | Temperature Limitation                                                                            |
|  | Humidity Limitation                                                                               |
|  | Atmospheric pressure limitation                                                                   |
|  | Class II                                                                                          |
|  | Waste must be controlled according to local regulations.                                          |
|  | Keep away from rain                                                                               |

|                                                                                     |                                                                                                                                      |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
|    | Non-Sterile                                                                                                                          |
|    | Do not use if package is damaged.                                                                                                    |
| <b>IP21</b>                                                                         | Ingress Protection<br>IP2X: Protected against solid objects over 12.5mm<br>IPX1: Protected against vertically falling drops of water |
|    | Universal serial bus (USB) port                                                                                                      |
|    | “ON” (power)                                                                                                                         |
|    | “OFF” (power)                                                                                                                        |
|    | DC Jack connector                                                                                                                    |
| <b>Rating</b>                                                                       |  12V 3A max                                         |
|  | Storage and Transportability Temperature Limitation                                                                                  |
|  | Storage and Transportability Humidity Limitation                                                                                     |
|  | Fragile, handle with care                                                                                                            |
|  | This Side UP Symbol<br>During of shipping or delivery, the carton should be kept upright.                                            |
|  | Stacked Boxes Maximum carton stack is 4.                                                                                             |
|  | Device Transport position                                                                                                            |
| <b>Rx ONLY</b>                                                                      | A medical prescription.                                                                                                              |

## 6.2. Operational Training

The device should be operated by a well-trained user. An on-site training for the Antlia Pro device is provided by a Customer Service Engineer or a 3rd party service engineer assigned by iCaltech during system installation. This training covers basic system operation and safety handling.

The following training should be provided to users:

1. Device Handling: Proper techniques for handling and transporting the device.
2. FOT Tube Positioning: Adjusting the height of the FOT Tube as per the patient height.
3. Usage of Single-use Accessories: Correct use of nose clips and iCaltech filters.
4. Device Verification Procedure: Steps involved in verifying the device accuracy.
5. Patient Positioning and Device Operation
6. Connection of accessories while performing test and device verification
7. Software Operation: For the Device Verification, Patient Testing and Report Generation

Apart from the above, iCaltech will provide necessary documentation and training to service personnel for assisting in parts repair.

## 6.3. Operating Procedure

### 6.3.1. User Login and Device Connection

- First, turn on the power supply to the device. Then, connect the device to the PC in which Antlia pro application installed, by using Type B USB cable. Then open the Antlia pro application.
- Note: The Antlia Pro application is provided by iCaltech and will be installed on the user's desktop or laptop by our representative during the device installation, following the instructions outlined in Chapter 7 of this manual.



Figure 12: Device in Ready to use mode

- Enter valid Username and Password with valid credentials and click "Login."



Figure 13: Antlia Pro Application Login Page

- The Home screen will be displayed on your PC as shown in figure 14.
- Click on the “Connect” icon, located in the top right corner of the home screen.




Figure 14: Antlia Pro Home Screen

- When the device is connected to the application, the status will change to "Connected", as shown in figure 15.



Figure 15: Device Connection Status on Home Screen

### 6.3.2. Device Verification Test

Each Antlia Pro device undergoes rigorous quality control to ensure maximum measurement accuracy. The sensors and electronics within the Antlia Pro are designed for excellent long-term stability. However, daily field verification is crucial to maintain system accuracy and promptly detect any potential changes or damage that may have inadvertently occurred.

This verification can be conducted by comparing the device's measurement to a known impedance value. Each Antlia Pro device is supplied with a specific

reference impedance device for this verification. The measured impedance value should be within a 10% tolerance of the known reference value.

### Device Verification Procedure:

- Connect the reference impedance to the FOT tube as shown in figure 16.



Figure 16 Reference Impedance connected to FOT tube

- Arrange the device and log in to the Antlia Pro application and connect the device with application as illustrated in section 6.3.1.
- Click on “Verification” button followed by the “Start” button. The device will initiate the verification process. The impedance value will begin displaying in the box on the left-hand corner under “Z(kPa/L/Sec)” after a few seconds as shown in the figure 17.



Figure 17: Device verification - Impedance Value Check

- If the value appears in “BLACK” color, it indicates that the device accuracy is within the acceptable accuracy limit. The user can then proceed with patient

testing.

- If the value appears in “**RED**” color, it indicates that the device accuracy is outside the acceptable limit. In this case, don’t proceed with the patient testing and contact the iCaltech Service team for assistance.

### 6.3.3. Performing Patient Tests

During device verification testing, if the impedance values fall within the acceptable accuracy limits, proceed to patient testing. Patient testing comprises four major steps:

1. Patient Information Entry
2. Patient Positioning
3. Performing the Test
4. Results and Reporting

#### 6.3.3.1. Patient Information Entry

- Go to Home screen on the application page → Click 'Patient Information' → 'Add Patient' (as shown in Figure 18).




Figure 18: Patient Information Page

- Enter all patient information in the given fields, as shown in the figure 19.



Figure 19: Patient Information Entry Page

- After entering the patient details, click the "Save and Conduct Test" button to save the patient details and begin the test immediately.
- If you wish to add details for multiple patients and conduct their tests individually later, click the "Save and Add More" button to continue adding patient records. Refer to figure 19 for a visual reference.

### 6.3.3.2. Patient Positioning

- Proper patient posture is crucial for accurate measurements.
- Place the device on a table as shown in the figure 20.



Figure 20: Positioning of Antlia Pro for Patient Test

- To perform the test correctly, position a chair in front of the table where the device is installed.
- Seat the patient in front of the device, ensuring they maintain an upright posture.
- Connect the iCaltech filter to the FOT tube of the device as illustrated in figure 21.



Figure 21: iCaltech Filter Connection

- Adjust the position of FOT tube to the patient mouth level by using the adjustable arm of the device.
- The technician must exercise caution with the source tube while adjusting the stand to prevent damage.
- Instruct patients to wear nose clips to prevent air leakage.
- During the test run, it is essential for the patient to hold their cheeks with their hands as shown in the figure 22 to enhance the reliability of the measurement.



Figure 22: Patient Positioning

- The patient should keep their chin up and breathe normally throughout the test.
- Ensure that there are no air leaks around the lips and iCaltech filter.

#### 6.3.3.3. Performing the Test

The Antlia Pro software is designed to assess respiratory function in accordance with the 2019 ERS guideline for Oscillometry. Therefore, a complete Antlia Pro patient test requires at least three separate, valid measurements. The overall test results are calculated as the average of these three tests.

##### **Measurement Modes:**

**Mono Frequency:** The PFT test should be carried out for 12 seconds at a frequency range of 5 to 19 Hz.

**Multi Frequency Mode:** The PFT test should be carried out in mixed frequency mode with a frequency range of 5, 11, 13, 17, 19, 23, 29, 31, and 37 Hz.

**Coefficient of Variation (COV) Limits:** In Mixed mode, the COV value should be less than 10% for adults and less than 15% for pediatrics for all three attempts. If the COV value exceeds these limits, the patient should repeat the test.

##### **Test Stages:**

The App supports conducting tests in two stages:

1. **Pre-test:** The Pre-test is selected when the patient is undergoing the test for the first time.
2. **Post-test:** Following the Pre-test, and based on the results, the physician may administer a bronchodilator, and then the Post-test is conducted.

##### **Performing the Test:**

- After entering the patient information and position the patient for the testing, click on “Save and Conduct the Test”.
- If patient information already exists in the database, the user has to select that patient from the conduct test screen.
- On the Home screen, Click the "Conduct test" button. This will display a list of patients added in the database. (Refer figure 23 for the visual reference).




Figure 23: Patient Database on Conduct Test Page

- Select the desired patient from the list and click the "Conduct test" button located at the bottom of the page (as shown in Figure 24).




Figure 24: Selecting Patient for Testing

- By default, the app will display the Pre-test option with Multi-mode selected as shown in figure 25.



Figure 25: Conduct Test Page

- The technician can select the Pre-test or Post-test option and the Multi or Mono-mode option as per their requirements.
- Select the Frequency mode ("Multi" or "Mono") from the dropdown box as shown in the figure 26.




Figure 26: Frequency Mode Selection

- Click the "Start Test" button and wait for the "Start Breathing Test" instruction to appear on the screen.




Figure 27: Start Breathing Test Instruction

- Instruct the patient to breathe through their mouth while gently holding their cheeks together to ensure a proper seal around the iCaltech filter and prevent any air leakage.
- Whenever the patient's breathing synchronizes with the waveform, click the "Capture" button.



Figure 28: Volume and Mouth Pressure Measurement

- After test completion, click the "View Report" button to display the test report on the screen.



Figure 29: Test Report Page

- Click the "Save as PDF" button to save the patient report of the current test.
- By clicking the "Test again" button, perform 3 trials with a COV value less than 10% for adults or less than 15% for pediatrics.

### 6.3.3.4. Results and Reporting

**Real time graph:** When the patient is breathing, we can see volume and mouth pressure graphs in the App.




Figure 30: Real Time Graph of Volume and Mouth Pressure Measurement


After completing the pre-test, 'Save' button must be pressed. The report will then be stored on the PC in pdf format at user defined location. After conducting post-test for the same patient, the 'Save' button must be pressed again. A consolidated report of the pre-



**To save the report, follow the procedure below:**

- [illegible]

- Save the PDF file to the desired location on the PC.



- The user can take a printout of the patient report from the PC.
- Go to the location where the report was saved. Open the PDF file. The report will look as depicted below.

## Page 47 of 72




**I. Patient with normal condition:****Respiratory and Healthy breath**

Bangalore

|                     |                    |                  |                       |
|---------------------|--------------------|------------------|-----------------------|
| Patient Name : Test | Patient ID : 10030 | Gender: Female   | Smoking status: Never |
| D.O.B : 01-Jan-99   | Age : 25 years     | Height: 163.0 cm | Weight : 72.0 kg      |
| BMI : 27.1          | Conducted by: sri  | Doctor: Dr.Latha | Test Date : 18-Dec-24 |

**RESPIRATORY OSCILLOMETRY**

| Parameters           | Cut Off* | PRE             | * ΔI (pre) | POST | * ΔI (post) | %Change | Bronchodilator % Change |   |   |
|----------------------|----------|-----------------|------------|------|-------------|---------|-------------------------|---|---|
| Z5 (mult, O/L/s)     | 4        | 3.22 (COV: 1.7) | -20 (s)    | -    | -           | -       | RS                      | - | - |
| Z10 (mult, O/L/s)    | -        | 2.27            | -          | -    | -           | -       | RS-R10                  | - | - |
| R5 (mult, O/L/s)     | 4        | 3.14 (COV: 1.5) | -22 (s)    | -    | -           | -       | X5                      | - | - |
| R5Sep (mult, O/L/s)  | -        | 3.06            | -          | -    | -           | -       | AX                      | - | - |
| R5Exp (mult, O/L/s)  | -        | 3.69            | -          | -    | -           | -       | Fms                     | - | - |
| R10 (mult, O/L/s)    | 3        | 2.15            | -28 (s)    | -    | -           | -       |                         |   |   |
| RS-R10 (mult, O/L/s) | 1        | 0.99            | -1 (s)     | -    | -           | -       |                         |   |   |
| X5 (mult, O/L/s)     | -1       | -0.71           | 29 (s)     | -    | -           | -       |                         |   |   |
| X5Sep (mult, O/L/s)  | -        | -1.27           | -          | -    | -           | -       |                         |   |   |
| X5Exp (mult, O/L/s)  | -        | -0.09           | -          | -    | -           | -       |                         |   |   |
| AX (mult, O/L/s)     | 4        | 3.14            | -          | -    | -           | -       |                         |   |   |
| Fms (Hz)             | <12      | 12.51           | -          | -    | -           | -       |                         |   |   |

**R and X graph:**

Printed on: 18/12/2024 15:47:56 with iCaltech Innovations Antlia Pro device \*Sah's et al ET AL research paper 2023

Page 1 of 3

## Respiratory and Healthy breath

Bangalore

|                                            |                    |                   |                        |
|--------------------------------------------|--------------------|-------------------|------------------------|
| Patient Name : Srinithi One<br>ninety four | Patient ID : 10030 | Gender : Female   | Smoking status : Never |
| D.O.B : 01-Jan-99                          | Age : 25 years     | Height : 163.0 cm | Weight : 72.0 kg       |
| BMI : 27.1                                 | Conducted by: sri  | Doctor : Dr.Latha | Test Date : 18-Dec-24  |

| Breathing Parameter  | PRE | POST |
|----------------------|-----|------|
| Resp. Rate/min       | 26  | -    |
| V <sub>T</sub> (ml)  | 355 | -    |
| T <sub>exp</sub> (s) | 1   | -    |
| T <sub>exp</sub> (s) | 1   | -    |


| Pre-Test Table<br>Parameters | Cut Off | Act 1 | Act 2 | Act 3 | Avg             | *AI     |
|------------------------------|---------|-------|-------|-------|-----------------|---------|
| R5 (cmH2O/L/s)               | 4       | 3.09  | 3.19  | 3.15  | 3.14(COV: 1.5)  | -22 (n) |
| R11 (cmH2O/L/s)              | 4       | 2.24  | 2.28  | 2.26  | 2.26            | -43 (n) |
| R19 (cmH2O/L/s)              | 3       | 2.18  | 2.15  | 2.13  | 2.15            | -28 (n) |
| R5-19 (cmH2O/L/s)            | 1       | 0.91  | 1.04  | 1.02  | 0.99            | -13 (n) |
| X5 (cmH2O/L/s)               | -1      | -0.65 | -0.69 | -0.6  | -0.71           | 29 (p)  |
| X11 (cmH2O/L/s)              | -       | -0.23 | -0.3  | -0.27 | -0.27           | -       |
| X19 (cmH2O/L/s)              | -       | 0.76  | 0.69  | 0.73  | 0.73            | -       |
| Z5 (cmH2O/L/s)               | 4       | 3.10  | 3.26  | 3.25  | 3.22(COV: 1.71) | -20 (n) |
| Z19 (cmH2O/L/s)              | -       | 2.31  | 2.28  | 2.25  | 2.27            | -       |
| Delta X5 (cmH2O/L/s)         | -       | 0.93  | 0.98  | 1.27  | 1.06            | -       |
| AX (cmH2O/L)                 | 4       | 2.70  | 3.05  | 3.4   | 3.14            | -       |
| CO @ 5Hz                     | -       | 0.03  | 0.02  | 0.03  | -               | -       |
| CO @ 19Hz                    | -       | 0.06  | 0.06  | 0.06  | -               | -       |
| Fres (Hz)                    | <12     | 12.2  | 12.73 | 12.61 | 12.51           | -       |
| X5insp (cmH2O/L/s)           | -       | -1.13 | -1.3  | -1.37 | -1.27           | -       |
| X5Exp (cmH2O/L/s)            | -       | -0.02 | 0.09  | -0.35 | -0.09           | -       |
| R5insp (cmH2O/L/s)           | -       | 3.07  | 3.16  | 2.94  | 3.06            | -       |
| R5Exp (cmH2O/L/s)            | -       | 3.41  | 3.67  | 3.96  | 3.68            | -       |
| X5EndInsp (cmH2O/L/s)        | -       | -1.2  | -1.56 | -1.42 | -1.39           | -       |
| X5EndExp (cmH2O/L/s)         | -       | 0.02  | -0.03 | -0.94 | -0.98           | -       |
| R5EndInsp (cmH2O/L/s)        | -       | 2.89  | 2.83  | 2.67  | 2.8             | -       |
| R5EndExp (cmH2O/L/s)         | -       | 3.39  | 3.6   | 3.27  | 3.42            | -       |
| Test Time                    | -       | 15:44 | 15:45 | 15:47 | -               | -       |

## Respiratory and Healthy breath

Bangalore

|                             |             |                    |                   |                        |
|-----------------------------|-------------|--------------------|-------------------|------------------------|
| Patient Name : Srinithi One |             | Patient ID : 10030 | Gender : Female   | Smoking status : Never |
| nienty four                 |             |                    |                   |                        |
| D.O.B                       | : 01-Jan-99 | Age                | : 25 years        | Height : 163.0 cm      |
| Weight                      |             |                    |                   | : 72.0 kg              |
| BMI                         | : 27.1      | Conducted by: sri  | Doctor : Dr.Latha | Test Date : 18-Dec-24  |

### Device Interpretation:



\*Normal

\* (") refers to Increment \* (n) refers to Normal \* (n) refers to Decrement \* AI refers to Anomaly Index \* Interpretation is done w.r.t. cut off values

### Doctor's Inference:

Test conducted by

Doctor's Signature

## II. Patient with abnormal condition:

### Respiratory and Healthy breath

Dr.Vinayak

MBBS


Kolkata

|                         |                      |                    |                        |
|-------------------------|----------------------|--------------------|------------------------|
| Patient Name : Samyukta | Patient ID : 700061  | Gender : Male      | Smoking status : Never |
| D.O.B : 01-Jan-65       | Age : 58 years       | Height : 162.0 cm  | Weight : 70.0 kg       |
| BMI : 26.67             | Conducted by : Durga | Doctor : Dr.Naveen | Test Date : 12-Dec-23  |

#### RESPIRATORY OSCILLOMETRY

| Parameters                                  | Cut Off | PRE              | * AI (pre) | POST            | * AI (post) | %Change | Bronchodilator % Change |        |      |
|---------------------------------------------|---------|------------------|------------|-----------------|-------------|---------|-------------------------|--------|------|
| Z5 (cmH <sub>2</sub> O/L/s)                 | 4       | 7.75 (COV: 10.0) | 94 (%)     | 6.23 (COV: 5.2) | 56 (%)      | -19.7   | RS                      | -15.5  | Decr |
| Z19 (cmH <sub>2</sub> O/L/s)                | 4       | 3.11             | -22 (n)    | 2.91            | -27 (n)     | -6.21   | RS-RIP                  | -34.4  |      |
| RS (cmH <sub>2</sub> O/L/s)                 | 4       | 6.77 (COV: 7.3)  | 69 (%)     | 5.72 (COV: 5.4) | 43 (%)      | -15.3   | AS                      | 34.32  | Incr |
| RS <sub>Stop</sub> (cmH <sub>2</sub> O/L/s) | -       | 5.74             | -          | 5.03            | -           | -12.3   | AX                      | -32.00 |      |
| RS <sub>Exp</sub> (cmH <sub>2</sub> O/L/s)  | -       | 8.76             | -          | 7.52            | -           | -14.15  | Fin                     | -2.02  |      |
| R19 (cmH <sub>2</sub> O/L/s)                | 3       | 3                | 0 (%)      | 2.87            | -4 (n)      | -4.33   |                         |        |      |
| RS-R19 (cmH <sub>2</sub> O/L/s)             | 1       | 3.76             | 276 (%)    | 2.84            | 164 (%)     | -24.4   |                         |        |      |
| Z5 (cmH <sub>2</sub> O/L/s)                 | -1      | -3.76            | -276 (v)   | -2.46           | -146 (v)    | 34.52   |                         |        |      |
| XStop (cmH <sub>2</sub> O/L/s)              | -       | -4.0029          | -          | -3.546          | -           | 16.41   |                         |        |      |
| XStop (cmH <sub>2</sub> O/L/s)              | -       | -3.9414          | -          | -1.5187         | -           | 61.46   |                         |        |      |
| AX (cmH <sub>2</sub> O/L/s)                 | 4       | 32.73            | -          | 22.23           | -           | -32.09  |                         |        |      |
| Fin (Pa)                                    | <12     | 24.4             | -          | 23.9            | -           | -2.02   |                         |        |      |

R and X graph:



Printed on: 12/12/2023 16:19:29 with iCaltech Innovations Antlia Pro device/Sahu's et al ET AL, research paper 2023, Only for adults > 18 years Page 1 of 4

## Respiratory and Healthy breath

Kolkata

|                         |                     |                    |                        |
|-------------------------|---------------------|--------------------|------------------------|
| Patient Name : Samyukta | Patient ID : 700061 | Gender : Male      | Smoking status : Never |
| D.O.B : 01-Jan-65       | Age : 58 years      | Height : 162.0 cm  | Weight : 70.0 kg       |
| BMI : 26.67             | Conducted by: Durga | Doctor : Dr.Naveen | Test Date : 12-Dec-23  |

| Breathing Parameter   | PNE   | POST  |
|-----------------------|-------|-------|
| Resp. Rate/min        | 22.67 | 24.33 |
| V <sub>T</sub> (ml)   | 624   | 547   |
| T <sub>insp</sub> (s) | 1.18  | 1.15  |
| T <sub>exp</sub> (s)  | 1.44  | 1.31  |

Pre-Test Table

| Parameters            | Cut Off | Act 1 | Act 2 | Act 3 | Avg              | 'AI      |
|-----------------------|---------|-------|-------|-------|------------------|----------|
| R5(cmH2O/L/s)         | 4       | 7.17  | 6.91  | 6.22  | 6.77(COV: 7.3)   | 69 (%)   |
| R11(cmH2O/L/s)        | 4       | 4.67  | 4.99  | 5.01  | 4.89             | 22 (%)   |
| R19(cmH2O/L/s)        | 3       | 2.88  | 3.04  | 3.09  | 3                | 0 (%)    |
| R5-19(cmH2O/L/s)      | 1       | 4.29  | 3.87  | 3.13  | 3.76             | 276 (%)  |
| X5(cmH2O/L/s)         | -1      | -4.68 | -3.88 | -2.93 | -3.76            | -276 (v) |
| Z5(cmH2O/L/s)         | 4       | 8.56  | 7.83  | 6.87  | 7.75(COV: 10.88) | 94 (%)   |
| Z19(cmH2O/L/s)        | 4       | 3.09  | 3.13  | 3.1   | 3.11             | -22 (v)  |
| AX (cmH2O/L)          | 4       | 43.96 | 32.81 | 21.42 | 32.73            | -        |
| CO @ 9Hz              | -       | 0.96  | 0.88  | 0.88  | -                | -        |
| CO @ 20Hz             | -       | 0.96  | 0.96  | 0.96  | -                | -        |
| Freq (Hz)             | <12     | 26.22 | 24.93 | 22.04 | 24.4             | -        |
| Test Time             | -       | 15.13 | 15.14 | 15.15 | -                | -        |
| X5insp (cmH2O/L/s)    | -       | -4.44 | -3.67 | -3.9  | -4               | -        |
| X5Exp (cmH2O/L/s)     | -       | -6.78 | -3.54 | -1.5  | -3.94            | -        |
| R5insp (cmH2O/L/s)    | -       | 8.08  | 5.5   | 5.63  | 5.74             | -        |
| R5Exp (cmH2O/L/s)     | -       | 9.54  | 8.31  | 8.42  | 8.76             | -        |
| X5EndInsp (cmH2O/L/s) | -       | -3.95 | -3.31 | -3.19 | -3.48            | -        |
| X5EndExp (cmH2O/L/s)  | -       | -9.53 | -5.6  | -3.04 | -6.06            | -        |
| R5EndInsp (cmH2O/L/s) | -       | 5.25  | 5.01  | 4.93  | 5.07             | -        |
| R5EndExp (cmH2O/L/s)  | -       | 8.37  | 9     | 9.11  | 8.83             | -        |


## Respiratory and Healthy breath

Kolkata

Patient Name : Samyukta Patient ID : 700061 Gender : Male Smoking status : Never  
D.O.B : 01-Jan-65 Age : 58 years Height : 162.0 cm Weight : 70.0 kg  
BMI : 26.67 Conducted by: Durga Doctor : Dr.Naveen Test Date : 12-Dec-23

| Post-Test-Table<br>Parameters | Cut Off | Act 1 | Act 2 | Act 3 | Avg             | *AI      | %chg   |
|-------------------------------|---------|-------|-------|-------|-----------------|----------|--------|
| R5(cmH2O/L/s)                 | 4       | 6.07  | 5.59  | 5.5   | 5.72(COV: 5.37) | 43 (*)   | -15.5  |
| R15(cmH2O/L/s)                | 4       | 4.83  | 4.69  | 4.37  | 4.63            | 16 (*)   | -5.34  |
| R19(cmH2O/L/s)                | 3       | 2.91  | 2.84  | 2.87  | 2.87            | -4 (n)   | -4.33  |
| R5-19(cmH2O/L/s)              | 1       | 3.16  | 2.74  | 2.83  | 2.84            | 184 (*)  | -24.4  |
| X5(cmH2O/L/s)                 | -1      | -2.56 | -2.5  | -2.33 | -2.46           | -146 (v) | 34.52  |
| Z5(cmH2O/L/s)                 | 4       | 6.59  | 6.12  | 5.97  | 6.23(COV: 5.18) | 56 (*)   | -19.7  |
| Z19(cmH2O/L/s)                | 4       | 2.96  | 2.9   | 2.88  | 2.91            | -27 (n)  | -6.21  |
| AX (cmH2O/L)                  | 4       | 23.61 | 24.51 | 18.56 | 22.23           | -        | -32.09 |
| CO @ 5Hz                      | -       | 0.92  | 0.86  | 0.85  | -               | -        | -      |
| CO @ 20Hz                     | -       | 0.97  | 0.97  | 0.97  | -               | -        | -      |
| Fres (Hz)                     | <12     | 24.1  | 24.36 | 23.24 | 23.9            | -        | -2.02  |
| Test Time                     | -       | 15:39 | 15:40 | 15:41 | -               | -        | -      |
| X5Insp (cmH2O/L/s)            | -       | -3.35 | -3.38 | -3.3  | -3.35           | -        | 16.41  |
| X5Exp (cmH2O/L/s)             | -       | -2.17 | -1.59 | -0.6  | -1.52           | -        | 61.46  |
| R5Insp (cmH2O/L/s)            | -       | 5.37  | 5.11  | 4.61  | 5.03            | -        | -12.3  |
| R5Exp (cmH2O/L/s)             | -       | 8.05  | 7.56  | 6.94  | 7.52            | -        | -14.15 |
| X5EndInsp (cmH2O/L/s)         | -       | -3.1  | -3.03 | -2.99 | -3.04           | -        | 12.58  |
| X5EndExp (cmH2O/L/s)          | -       | -3.38 | -3    | -1.59 | -2.66           | -        | 56.13  |
| R5EndInsp (cmH2O/L/s)         | -       | 4.86  | 4.55  | 3.86  | 4.42            | -        | -12.69 |
| R5EndExp (cmH2O/L/s)          | -       | 8.5   | 8.02  | 7.56  | 8.02            | -        | -9.07  |

### Device Interpretation:



Larger Airway Obstruction : 100.11 %

Smaller Airway Obstruction : 376 %

Lung Elasticity : 376 %

\* (\*) refers to Increment \* (n) refers to Normal \* (v) refers to Decrement \* AI refers to Anomaly Index \* Interpretation is done w.r.t. cut off values

## Respiratory and Healthy breath

Kolkata

|                         |                     |                    |                        |
|-------------------------|---------------------|--------------------|------------------------|
| Patient Name : Samyukta | Patient ID : 700061 | Gender : Male      | Smoking status : Never |
| D.O.B : 01-Jan-65       | Age : 58 years      | Height : 162.0 cm  | Weight : 70.0 kg       |
| BMI : 26.67             | Conducted by: Durga | Doctor : Dr.Naveen | Test Date : 12-Dec-23  |

|                                 | Pre - Post % change | Absolute Change |
|---------------------------------|---------------------|-----------------|
| R5 (cmH <sub>2</sub> O/L/s)     | -15.5 %             | -1.1            |
| R19 (cmH <sub>2</sub> O/L/s)    | -4.33 %             | -0.1            |
| R5-R19 (cmH <sub>2</sub> O/L/s) | -24.4 %             | -1              |
| N5 (cmH <sub>2</sub> O/L/s)     | 34.52 %             | -1.3            |
| AN (cmH <sub>2</sub> O-L)       | -32.09 %            | -10.5           |

\* - Doctor's review comments required

Doctor's comments:

Test conducted by

Doctor's Signature



## 6.4. Care And Maintenance

### Accessories Disposal:

- iCaltech filters are for single-use only. Using them correctly helps prevent the transmission of diseases and infections between patients.
- Upon completion of a patient test, ensure proper disposal of single-use supplies, including the iCaltech filter, in a suitable waste collection container according to your institution's guidelines and regulations.

### Storage and Transportation:

- When not in use for extended periods, cover the device with the dust cover.
- For long-term storage or transportation, use the original Antlia Pro packaging.
- After storing or transporting in cold or warm environments, allow several minutes for the device to acclimatize to room temperature before use. This ensures accurate measurement results.

### Daily Maintenance:

- Perform daily device verification before each use (refer to section 6.2: Device Verification Test).

### Annual Maintenance:

- The Antlia Pro device requires annual service by iCaltech-authorized service providers.
- Contact iCaltech Technical Support or your local representative to schedule service.
- Interchangeable or detachable parts specified by iCaltech must be replaced by a service engineer.
- Component or part failures must be addressed according to the iCaltech Service Protocol.

### Contact Details:

- Toll-Free Number: 1800-599-3196

### Expected Service Life:

- The expected service life is five years from the date of purchase.

### Routine Cleaning:



- Proper use and disposal of disposable antibacterial filters and nose clips contribute to infection prevention by minimizing the risk of disease transmission between patients and reducing the need for extensive device cleaning.
- Regularly clean the exterior surfaces of the device with a mild detergent solution and a soft cloth.

**CAUTION**

Do not use harsh chemicals, abrasive cleaners, alcohol, or sterilization liquids on the device surface.

**Disposal of Device Post Service Life:**

- At the end of its service life, dispose of the device according to all local regulations for electrical and electronic waste.
- Improper disposal can harm the environment and human health.
- The crossed-out wheeled bin symbol indicates that the product must be taken to a designated collection point for proper disposal.

## 6.5. Troubleshooting

While using the device, users may encounter common operating problems or errors. Below is a list of the most frequent issues that may occur during normal use, along with their possible solutions. If you experience a problem not listed here, please contact the iCaltech Service Team for further assistance.

### List of Common Problems and Solutions:

| Sl. No. | Problem                                                    | Possible Cause                           | Solution                                                                                                                             |
|---------|------------------------------------------------------------|------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| 1       | Unit fails to Power ON; all power indicators remain unlit. | Power cable not connected                | Check if power cable is properly connected to the device and the power outlet                                                        |
|         |                                                            | No power at the outlet                   | Check if power is available at the outlet. Try plugging the device into a different outlet.                                          |
|         |                                                            | Power button is in the OFF condition     | Press the power ON button to turn on the device.                                                                                     |
| 2       | Device not connecting with APP                             | No power on outlet                       | Switch ON the device                                                                                                                 |
|         |                                                            | USB cable not properly connected         | Disconnect and reconnect the USB cable securely.                                                                                     |
|         |                                                            | PC not Compatible with the software      | Check if the COM port drivers are installed on your PC and confirm PC meets the specifications required for the software             |
| 3       | RED LED is glowing continuously                            | Software error                           | Check error messages on screen and proceed accordingly.                                                                              |
|         |                                                            | Internal error (rare)                    | Disconnect and connect power supply cable and Re-start the system                                                                    |
| 4       | Mouth pressure and Flow rate graph is not displaying       | Incorrect patient position               | Check patient is in correct posture, wearing nose clip and holding their cheeks. Check if the iCaltech filter is connected properly. |
|         |                                                            | Defective Flow sensor or Pressure sensor | Replace Flow sensor or pressure sensor, if defective.                                                                                |

| Sl. No. | Problem                                                                                                                         | Possible Cause                                       | Solution                                                                                                              |
|---------|---------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|
|         |                                                                                                                                 | Breathing circuit defective or incorrectly assembled | Verify Breathing circuit is defective and properly assembled. Check error messages on screen and proceed accordingly. |
| 5       | Device impedance is exceeding 10% accuracy during verification test by using reference impedance                                | Improper connection of Reference impedance to device | Connect the known resistance to FOT tube without iCaltech filter                                                      |
|         |                                                                                                                                 | Incorrect Known resistance is used                   | Check Known resistance value and perform the test                                                                     |
| 6       | Application crashes while switching from the 'Start' to the 'Capture' button                                                    | Internal error                                       | Restart the application software                                                                                      |
| 7       | Unable to save the report after pre/post-test completion                                                                        | Internal error                                       | Restart the application software                                                                                      |
| 8       | When application is closed, the connection between device and the application is disconnected, but the device continues to run. | Internal error                                       | Restart the device                                                                                                    |

### List of Errors and Actions

| Sl. No. | Error Code | Error Type   | Error message Title | Error Message                        | Scenario                                                                                                                            | Actions                                  |
|---------|------------|--------------|---------------------|--------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------|
| 1       | NA         | Text message | NA                  | Please enter <missing field name>    | If any of the mandatory field is missing in edit/ new patient screen, hospital screen, Staff screen and user clicks on save button. | Enter all the mandatory fields with *    |
| 2       | NA         | Text message | NA                  | Phone number can contain only digits | If phone number is entered in alphabets, alphanumeric.                                                                              | On phone number field enter only Numbers |

| Sl. No. | Error Code         | Error Type   | Error message Title             | Error Message                                                          | Scenario                                                                                                   | Actions                                          |
|---------|--------------------|--------------|---------------------------------|------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|--------------------------------------------------|
|         |                    |              |                                 |                                                                        | (Hospital information screen)                                                                              |                                                  |
| 3       | NA                 | Text message | NA                              | First name Must be contained only alphabets with maximum 20 characters | If First name field is invalid. (Add new staff screen)                                                     | Enter only Alphabet's with Max.20 characters     |
| 4       | NA                 | Text message | NA                              | Search not found                                                       | When user search not found in patient or staff screens. (Patient list screen)                              | Enter valid patient name or ID                   |
| 5       | NA                 | Text message | NA                              | Please enter valid username and password                               | If both or any one of username and Password fields are invalid. (Login screen)                             | Enter valid User name and password               |
| 6       | NA                 | Text message | NA                              | Please enter the username. Please enter the password.                  | If username field or password field is blank (Login screen)                                                | Do not keep username and password field as blank |
| 7       | NA                 | Text message | NA                              | New password and Confirm password must be same                         | If entered new password and confirm password are not same. (Change password screen. Admin Login)           | Enter valid password                             |
| 8       | AS_E<br>RR_0<br>05 | Information  | No reports                      | No patient reports are available                                       | If user clicks on view report in home screen, when test is not conducted for patient. (View report screen) | Conduct test and click on view report screen     |
| 9       | AS_E<br>RR_0<br>06 | Information  | No reports for selected patient | "No reports are available for selected patient"                        | If no reports are available for the selected patient. (View report screen)                                 | Conduct the test and save the Report             |
| 10      | AS_E<br>RR_1<br>02 | Error        | Patient not selected            | Please select patient from the list                                    | If user clicks on Begin test button without selecting patient from                                         | Select patient and click on Begin test           |

| Sl. No. | Error Code         | Error Type | Error message Title                      | Error Message                                                                                                | Scenario                                                                                                         | Actions                                        |
|---------|--------------------|------------|------------------------------------------|--------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|------------------------------------------------|
|         |                    |            |                                          |                                                                                                              | the list. (Home screen)                                                                                          |                                                |
| 11      | AS_E<br>RR_1<br>08 | Error      | Problem connecting with device           | Please check USB cable connection and click on "Connect" button                                              | Without connecting device to PC. If user searches for particular device calibration values. (Calibration screen) | Connect the device with APP before calibration |
| 12      | AS_E<br>RR_1<br>10 | Warning    | Hospital Information delete confirmation | "Are you sure, you want to delete Hospital information?" with "Delete" and "Cancel" options                  | After selecting hospital, if user clicks on delete button. (Hospital information screen)                         | Confirm about deletion of the filed            |
| 13      | AS_E<br>RR_1<br>11 | Warning    | Staff delete confirmation                | "Are you sure want to delete staff <staff full name>" with "Delete" and "Cancel" options                     | If user clicks on delete button of particular staff information. (Staff list screen)                             | Confirm about deletion of the filed            |
| 14      | AS_E<br>RR_1<br>12 | Warning    | Cancel confirmation: Add staff           | All the entered staff information will be lost                                                               | If user clicks on cancel button after entering or editing the staff details                                      | Confirm about cancelation of the filed         |
| 15      | AS_E<br>RR_1<br>13 | Warning    | Cancel confirmation: Edit Staff          | All the entered staff information will be lost<br>"Are you sure want to cancel?" with "Yes" and "No" options | If user clicks on cancel button after entering or editing the staff information. (Edit staff screen)             | Confirm about cancelation of the filed         |
| 16      | AS_E<br>RR_1<br>14 | Warning    | Delete Confirmation: Patient Information | Are you sure want to delete patient <Patient full name> with "Delete" and                                    | If user clicks on delete button of particular patient information. (Add patient screen)                          | Confirm about deletion Of the filed            |

| Sl. No. | Error Code         | Error Type  | Error message Title              | Error Message                                                                                                | Scenario                                                                                     | Actions                                          |
|---------|--------------------|-------------|----------------------------------|--------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|--------------------------------------------------|
|         |                    |             | on                               | "Cancel" options                                                                                             |                                                                                              |                                                  |
| 17      | AS_E<br>RR_1<br>15 | Warning     | Cancel confirmation: Add Patient | All the entered patient information will be lost. "Are you sure want to cancel?" with "Yes" and "No" options | If user clicks on cancel button after entering the patient information. (Add patient screen) | Confirm about cancelation of the filed           |
| 18      | AS_E<br>RR_1<br>16 | Warning     | Logout Confirmation              | Are you sure want to log out? with "Ok" and "Cancel" buttons                                                 | If user clicks on Log-out button (Home screen)                                               | Confirm about Logout                             |
| 19      | AS_E<br>RR_1<br>17 | Warning     | Exit Confirmation                | Are you sure want to exit from the Antlia Pro App? With 'Ok' and 'Cancel' buttons                            | If user close (X) the application (Home Screen)                                              | Confirm about closing the application            |
| 20      | AS_E<br>RR_1<br>18 | Information | Duplicate staff Information      | Staff information is already existing                                                                        | For duplicate staff creation. (Add new staff screen)                                         | Do not enter already Entered staff information.  |
| 21      | AS_E<br>RR_1<br>19 | Information | Duplicate Patient Information    | Patient information is already existing                                                                      | For duplicate patient creation. (Add patient screen)                                         | Do not enter already entered patient information |
| 22      | AS_E<br>RR_1<br>20 | Information | Maximum trials reached           | Test has reached to maximum trials, please save the report before continuing with new test.                  | While performing the 10th action in pre/post mode. (View report screen)                      | After 10th Action save the Report                |

| Sl. No. | Error Code         | Error Type  | Error message Title    | Error Message                                                          | Scenario                                                                                                     | Actions                                                    |
|---------|--------------------|-------------|------------------------|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|------------------------------------------------------------|
| 23      | AS_E<br>RR_1<br>21 | Information | Maximum trials reached | Test has reached to maximum trials, Do you want to replace the values? | While performing the 11th action in pre/post mode. (View report screen)                                      | After 10th Action save the Report                          |
| 24      | AD_E<br>RR_2<br>02 | Error       | Device Disconnection   | Please connect the device to conduct the test                          | In between if user disconnects USB cable from the PC. (Home screen)                                          | Do not disconnect the USB cable, while performing test     |
| 25      | AD_E<br>RR_2<br>03 | Error       | Device disconnected    | Please connect the device to conduct the test                          | During test, if device got reset due communication slow down. (Conduct test screen)                          | Disconnect the USB cable and reconnect                     |
| 26      | AD_E<br>RR_2<br>04 | Error       | Device Disconnection   | Please connect the device to conduct the test                          | If user power off the device when USB is connected with PC. (Home screen)                                    | Do not power off the device while performing test          |
| 27      | AD_E<br>RR_2<br>08 | Error       | Device disconnected    | Please connect the device to conduct the test                          | While communicating with Device in Ideal mode                                                                | Disconnect the USB cable an reconnect and perform the test |
| 28      | AD_E<br>RR_2<br>09 | Error       | Device Accuracy Low    | Requires re- calibration                                               | Device impedance above 10% accuracy using known impedance value. (Verification screen – during verification) | Perform Device calibration using reference impedance       |
| 29      | AD_E<br>RR_2<br>11 | Warning     | Timeout                | Capture timeout. Please start the test again when ready                | Start mode running more than 2mins. (Conduct test screen)                                                    | Stop the test and perform the test again                   |

## 7. Application Software Design

### 7.1. Device connection with APP

#### 7.1.1. Software Installation

Antlia Pro Desktop application will be installed by iCaltech team. Install the application software according to user's PC's configuration.

- Note: System requirements: Operating System: Windows 10.0 or later & System Architecture: 64 bits

##### 7.1.1.1. Installation steps

###### Step 1: Software File

- Insert the Pen drive/CD that contains the Antlia Pro software installation .exe file.
- The contents will be as shown in Figure 33.




Figure 33: Antlia Pro Installation File

###### Step 2: Check PC's Specifications

- Right-click on "This PC" icon in your PC as shown in Figure 34 and click on "Properties."






Figure 34: PC's Properties

- This will open the system properties window and display the system information.
- In this screen, verify the operating system type (64-bit or not) as shown in Figure 35.




Figure 35: System Type - Verification

### Step 3: Start the Installation

- Double-click on the Antlia Pro software installation .exe file.

- A dialog box will appear requesting permission to install the app. Click "Yes."

#### Step 4: Select the Location of Installation Files

- Select the location where you want to install the application.
- The default location is "Program Files," but you can click on "Browse" to choose a different location. Refer Figure 36 for the location of installation folder.
- After selecting the location, click "Next."




Figure 36: Installation Folder Location in PC

#### Step 5: Desktop Shortcut Icon

- Check the checkbox if you want to create a desktop shortcut for easy access to the application. Refer figure 37 for the details.
- Otherwise, leave the checkbox unchecked → Click "Next."
- Note: Even now, if you want to change the installation location, click on "<Back>" and repeat step 4.




Figure 37: Selecting for Desktop Shortcut Icon

### Step 7: Install the Application

- Click "Next" to proceed with the installation.
- Wait for the installer to complete the process. The dialog box will disappear when the installation is finished. (Image need to be changed)




Figure 38: Installation Process - Completed

### Step 8: Launch the Application

- If a desktop shortcut was created: Double-click on the Antlia Pro desktop icon.



Figure 39: Antlia Pro - Desktop Icon

- If no desktop shortcut was created: Navigate to the installation directory → Double-click on "Antlia Pro-Vx" to launch the application.
- Upon launching the app, a splash screen will be displayed and after that the main application screen will appear as shown in figure 40.



Figure 40: Antlia pro - Main application screen

### 7.1.2. Software Upgradation

- If the user wants to upgrade from the old version to new version, install the new Antlia Pro application on your PC following steps outlined in 7.1.2.
- After that, navigate to the folder where old version is installed (eg., typically C:\Program Files\AntliaPro-Va.bb.ccd-Xee\App).
- Copy the database file named “antlia\_local.db” from the old application folder to the new application folder.

- Launch the newly installed Antlia pro application and proceed with testing. Refer to figure 41 for the details of data file location if needed.




Figure 41: Data base file location in program files

### 7.1.3. Firmware Upgradation

- Login to the admin window and navigate to the Firmware Upgrade screen.



Figure 42: Firmware Upgradation Page

- Observe the device's LED status.
- Verify if the device is already in firmware upgrade mode (Indicated by Blue and Red light glowing on the device).
- If the device is not in firmware upgrade mode (green LED is glowing), click "NO" option in the firmware upgradation page.
- Once the "No" option is selected, the "Enter Firmware Upgrade Mode" option will appear on the page as shown in figure 43.




Figure 43: Firmware upgrade mode option

- Click on "Enter Firmware Upgrade Mode." The device will then enter firmware upgrade mode, and the Blue and Red LEDs will start glowing.
- Now select “Yes” for the prompt “Is the device in firmware upgrade mode?”.
- The “Choose Firmware Hex file” option will appear on the page as shown in figure 44.
- Click on it, browse to the location of the Firmware Hex File, and select it.
- After that, the "Start Upgrade” option will appear on page as shown in figure 45. Click on it.
- Wait until a “success” pop-up message is displayed as shown in figure 45, indicating that the firmware upgrade has been completed successfully.



Figure 44: Choose Firmware HEX File




Figure 45: Start Upgrade option on firmware upgrade page

## **8. Exit from Application**

### **8.1. Software Closure**

To close the Antlia Pro software, click the close (X) button located in the top right corner of the Antlia Pro screen.

### **8.2. Antlia Pro Device Power Off**

Once usage of device finished, press and hold the Power ON/OFF button on the back of the device for approximately 3 seconds. Wait for all LED indicators on the device to turn off completely. Switch off the main power supply by turning the Main's switch on the Mains Board to the "OFF" position.



## Warranty Certificate

A sample Warranty card is shown below. The actual warranty card will be provided along with the device

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | <b>ICALTECH INNOVATIONS PRIVATE LIMITED</b><br>2nd Floor, Sy no 80, Municipal No 4 & 4A (Katha no 351),<br>Muneshwara Layout, Service Road, Ring Road,<br>Laggere Bridge, Laggere, Bengaluru, Karnataka, 560058, India.<br>Ph. # +91-94800 17633; www.icaltech.com<br>CIN: U74999KA2016PTC094056 GST: 29AAGCC5274M1Z3 |
| <p style="text-align: center;"><b><u>Product Warranty Certificate</u></b></p> <p><b>Customer Name:</b> &lt;Customer/Hospital/Clinic Name&gt;<br/><b>Kind Attn:</b></p> <p>We thank you for your recent purchase of our product detailed here under:</p> <p><b>Product Name:</b> Respiratory Impedance Monitoring Device ("<b>&lt;Product Name&gt;</b>")<br/><b>Product Sale Details:</b></p> <p><b>LOT No.:</b> &lt;Product Lot Number&gt; ; <b>Serial No.:</b> &lt;Product Serial Number&gt;</p> <p><b>Invoice No.:</b> &lt;Invoice Number&gt;</p> <p><b>Warranty Duration:</b> &lt;MMM YYYY&gt; to &lt;MMM YYYY&gt;<br/><b>AMC/CMC if any</b> _____</p> <p>We are glad to provide you this Warranty Certificate for your "<b>&lt;Product Name&gt;</b>" purchase. "<b>&lt;Product Name&gt;</b>" has passed Rigorous Tests and has been confirmed to be of Excellent Quality and Construction.</p> <p>If the "<b>Customer</b>" has any issue with the purchase during the Warranty Period, they can approach "<b>iCaltech</b>" (Details provided at end of this Warranty Certificate). "<b>iCaltech</b>" will initiate necessary actions within 3 Business Days to remedy the situation. Actions may include analyzing / rectifying the issue over Phone / Emails, in person visit to "<b>Customer</b>" place, or finally replacing the product itself. "<b>iCaltech</b>" will try its best to minimize the lead time. However, all actions to rectify any issue with the Product will be at the discretion of "<b>iCaltech</b>". The period of Warranty will be restricted to the Dates mentioned above, even if the product has been replaced in between the Warranty Period.</p> <p><b><u>This Warranty does not cover damage or defects caused by:</u></b></p> <ul style="list-style-type: none"><li>• Improper usage (failure to follow guidelines/ Instructions on using the product)</li><li>• Negligence, accident, damage due to inefficiencies in packaging during shifting by the "<b>Customer</b>".</li><li>• Any attempt to OPEN the device, customize the product with Modifications, use of unauthorized parts, change of Design, connections to any External products without "<b>iCaltech</b>" authorization</li><li>• Use third party filters / consumables without authorized by "<b>iCaltech</b>".</li></ul> <p><b><i>iCaltech Innovations Private Limited</i></b> has been referred to as "<b>iCaltech</b>" in this document.<br/><b><i>Manipal Health Enterprises Pvt Ltd</i></b> has been referred to as "<b>Customer</b>" in this document.</p> <p><b><u>Business Days Definition:</u></b><br/>Monday to Friday: 9 am to 6 pm<br/>(<b>Note:</b> Saturday / Sunday / Government Holidays will not be considered Business Days)</p> |                                                                                                                                                                                                                                                                                                                       |